

VAC-Via Vacuum Assisted Closure (VAC) Therapy

What is V.A.C.® Therapy?

VAC Therapy helps a wound heal. The VAC Via Therapy unit is a system with a soft vacuum. The unit is battery operated. It:

- Helps close the wound.
- Removes fluids from the wound.
- Limits the need for daily dressing changes.

The system has five parts:

1. A therapy unit or machine.
2. Foam dressings. The dressing goes on or inside the wound.
3. Clear adhesive drape. The drape covers the dressing.
4. Plastic tubing. One end of the tubing connects to the dressing. The other end connects to the canister that fits into the unit.
5. Canister.

When the dressing is first put on, some children feel a mild pulling. This is normal. After a few minutes, this goes away.



What do I need to know about using the VAC-Via unit at home?

- Check the dressing. Be sure the foam dressing stays shrunk down. It should look like a wrinkled raisin.
 - If it does not look like a wrinkled raisin, call the clinic. If it is after hours, call the doctor on call or go to the Emergency Department.
- Make sure the machine is running when the dressing is on.
 - Never turn it off for more than two hours.
 - Do not change the settings.

Safety

- Be careful that no one trips on the extra tubing.
- Be careful that the tubing does not catch on anything.

System alarm

The unit can alarm for many reasons. To solve alarms:

- Call The KCI Product Issue Response team at 1-800-275-4524.
- If KCI cannot solve the problem, call the clinic during the day time if this occurs. After hours, call the doctor on call or go to the Emergency Department.

Activity

How much your child can move around depends on their condition and the reason for the VACVia Therapy.

Your child will wear the VAC Via Therapy unit using the belt or strap provided with the therapy unit.

See the After Visit Summary (AVS) for your child's activity level.

What do I do when we no longer need the unit?

- These parts of the VAC Via Therapy unit should be recycled:
 - VAC Via.
 - Power cord.
 - AC adapter.

A specially designed, safe, postage paid shipping container has been provided.

- Put only the VAC Via Therapy unit into the shipping container. Seal container. Put it in the mail.

ALERT: Call your child's doctor, nurse, or clinic if you have any concerns or if your child has:

- A fever of 101.5° F (38.6°C).
- More pain at the wound.
- Bloody drainage in the VAC tubing.
- Pus or a bad smell from the wound.
- More redness or warmth around the wound.
- Special health care needs not covered by this information.

Also call if:

- The foam dressing does not stay shrunk down and it does not look like a wrinkled raisin.
- You are having problems with the VAC Via Therapy machine.
- The VAC Via Therapy has been turned off for more than two hours.

For other health and wellness information, check out this resource:

<https://kidshealth.org/ChildrensWi/en/parents>

This sheet was created to help you care for your child or family member. It does not take the place of medical care. Talk with your healthcare provider for diagnosis, treatment and follow-up.